



BabyPod™ II

INFANT TRANSPORT DEVICE

Transporting nature's precious cargo

Features

- Cushioned and insulated to protect baby during transport
- Lightweight unit for one-operator transport
- Adaptable to many environments



@operSurgical

Keeping You at the Forefront
of Women's Health Care™

BabyPod™ II is designed to safely convey a pediatric patient between departments within the hospital or between medical facilities.

Protects Baby

- Provides thermo-stable environment.
- Shock absorbent foam absorbs vibration during transport.
- Vacuum mattress conforms to baby's shape to limit unwanted movement.
- Viewing lid protects baby from loose moving objects while allowing visual monitoring.

Easy to Use

- Lightweight carbon fiber shell
- 20-pound unit suitable for one-operator transport
- Straps quickly to any stretcher or gurney

Flexible

- Non-electric heat source: TransWarmer® Infant Transport Mattress.
- Infusion pumps and/or IV's attach directly to the unit.
- BabyPod™ II may be hand-carried over rough terrain or when a stretcher is unavailable.

Multiple uses

- Evacuation
- Search and rescue
- Ground ambulance
- Between hospital transport
- Inter-department transport



TransWarmer® mattress keeps infants warm while in the BabyPod™ II or, as a stand alone heat source in the NICU, nursery, or anywhere newborns can experience cold stress.



Only nature builds a better nest



TransWarmer® Infant Transport Mattress

TransWarmer® Infant Transport Mattress provides a thermostable environment during transport for up to two hours. The TransWarmer Infant Transport Mattress is a gel-filled, disposable, thermal mattress designed for use in neonatal transports. TransWarmer's patented WarmGel® filler provides even heating to minimize the risk of hypothermia.

Features & Benefits

Minimizes the risk of cold stress

- Patented WarmGel® mattress
 - Nests infant, absorbing vibration during transport
- Consistent, uniform heat distribution
 - Prevents the complications of cold stress
- Heats in less than 60 seconds
 - Provides rapid thermal support
- Easy to activate
 - Simple, no inner bags to rupture
- Non-toxic, food-grade gel

BabyPod™ II

Transporting nature's precious cargo

Infant Transport Device

BabyPod™ II consists of a lightweight carbon fiber outer shell, lined with a shock absorbent foam inner layer, and has transparent shields for viewing the baby. It also contains a patient positioning vacuum mattress, stretcher fixing straps to secure the device for transport, and safety straps to secure the infant inside the device during transport.

Vital Statistics

Device Dimensions	38" x 17" x 12"
Device Weight (No accessories)	20lbs/9kg
Patient Weight Maximum	Up to 20lbs/9kg



Ordering Information

70002	BabyPod™ II Infant Transport Device (1 each)
20421	TransWarmer® Infant Transport Mattress (6 per box)

BabyPod™ II



CooperSurgical

*Keeping You at the Forefront
of Women's Health Care™*

Transwarmer® is a registered trademark of Prism Enterprises, L.P.

95 Corporate Drive, Trumbull, CT 06611 • 800.243.2974 • 203.601.5200

www.coopersurgical.com • ©2006 CooperSurgical, Inc.

Transwarmer and WarmGel are registered trademarks of Prism Enterprises, a CooperSurgical company.
Form # 80641 Rev. 08/06