

# Infla10<sup>®</sup> AX

## inflatable penile prosthesis with Anatomical Expansion

### Product Overview

Infla10<sup>®</sup> AX is indicated for patients suffering from chronic, organic erectile dysfunction (impotence) and are candidates for a penile prosthesis implantation.<sup>1</sup>

Infla10<sup>®</sup> AX is currently the only inflatable penile prosthesis on the market offering length expansion and hydrophilic coating.

Infla10<sup>®</sup> AX is designed for you and your patients and offers innovative features<sup>1</sup>:

- ♂ **Length Expansion**
- ♂ **HydroShield™ Coating**
- ♂ **RapidPump™**
- ♂ **AdaptiveReservoir™**
- ♂ **EasyClick™ Connectors**
- ♂ **ConnectSecure™ Extenders**



# Key Features & Benefits

Infla10® AX is designed with years of experience in urology and with one thing in mind; an innovative device that offers safety and effectiveness for you and your patients.

## Length Expansion

Offers the patients anatomical girth and length expansion addressing concerns regarding post-operative penile shortening.

## HydroShield™ Coating

Infla10® AX incorporates a hydrophilic coating on all external component surfaces. The hydrophilic coating offers you the freedom to choose the aqueous solution. Hydrophilic coating facilitates the rapid and strong absorption of the solution on the device and may promote easier device implantation.

## RapidPump™

Our innovative pump mechanism offers greater ease and speed for your patients when inflating and deflating the prosthesis.

## AdaptiveReservoir™

The reservoir will easily take the shape of the implantation site offering post-op comfort for your patients.

## EasyClick™ Connectors

All tubing connects easily without any additional assembly tools. Simply click to secure connections.

## ConnectSecure™ Rear Tip Extenders (RTE)

ConnectSecure™ RTE provides a secure and rigid connection.

Infla10® AX Infrapubic Approach	
Prosthesis Product Code	Prosthesis Lengths
INF1012AX-IP	12 cm
INF1015AX-IP	15 cm
INF1018AX-IP	18 cm
INF1020AX-IP	20 cm
INF1022AX-IP	22 cm
INF1024AX-IP	24 cm

Infla10® AX Penoscrotal Approach	
Prosthesis Product Code	Prosthesis Lengths
INF1012AX-PS	12 cm
INF1015AX-PS	15 cm
INF1018AX-PS	18 cm
INF1020AX-PS	20 cm
INF1022AX-PS	22 cm
INF1024AX-PS	24 cm

Infla10® AdaptiveReservoir™	
Reservoir Product Code	Volume
INF10RS-65	65 ml
INF10RS-70	70 ml
INF10RS-110	110 ml

Infla10® ConnectSecure™ Rear Tip Extenders	
0.5 cm / 1 cm / 1.5 cm / 2 cm / 3 cm / 4 cm / 5 cm / 6 cm	



Infla10® AdaptiveReservoir™  
70ml - 110ml



Infla10® AdaptiveReservoir™  
65ml

Caution: Federal Law (U.S.) restricts this device to sale by or on the order of a physician.

Prior to use, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions and potential adverse events.

References: 1 Rigicon® Infla10® AX three-piece penile prosthesis Instructions for Use.

Rigicon, INC.

2805 Veterans Memorial Hwy STE 5, Ronkonkoma,

NY 11779 United States of America

+1 (888) 202-9790 - www.rigicon.com