


LAVH → TLH

CerviGrip + ColpoWave

Rotating edge gives you the edge




Adjustable tip ✓
fits all uterine sizes

Cervical fixation for ✓
predictable ureteric
and uterine distances

Two colpotomizer ✓
sizes in one product


Raised Margins for ✓
colpotomy and oncology



CerviGrip™

Uterine Manipulator

Secure Cervical Fixation with Adjustable Uterine Tip



General Laparoscopy

Adjustable uterine tip to fit all uterine lengths ✓

Built-in tenaculum holder ✓

Tight cervical seal for hydrotubation ✓

Syringe port for hydrotubation ✓

TLH / LAVH

Adjustable tip removes need for various tips ✓

Secure cervical fixation ✓

Predictable ureteric and uterine distances with ColpoWave ✓

Facilitates uterine removal ✓




SurgiToolsMIS.com

ColpoWave™

Vaginal Delineator

For TLH and Sacrocolpopexy



Use as Standalone

Two colpotomizer sizes in one product ✓

Raised Margins for colpotomy and oncology ✓

Maintains pneumoperitoneum ✓

Sacrocolpopexy platform ✓

Use with CerviGrip


✓ Adjustable tip removes need for various tips

✓ Secure cervical fixation

✓ Predictable ureteric and uterine distances

✓ Facilitates uterine removal

 **SurgiToolsMIS.com**




LapNeed™

Laparoscopic suturing made easy!

Simpler Suturing

4 / 5 lap suturing steps

- Needle Introduction ✓
- Needle Handling ✓
- Needle Transfer ✓
- Needle Removal ✓
- Knot Tying 

Control & Versatility

Suitable for all forms of MIS


- ✓ General Surgery
- ✓ Gynaecology
- ✓ Urology
- ✓ Bariatric Surgery
- ✓ Paediatric Surgery

SurgiAssist™

Delivering a stable platform for laparoscopic surgery

Uterine Positioner

- Provides stable pelvic platform ✓
- No power or hydraulics ✓
- Easy and flexible ✓
- Free up assistant ✓



Camera Holder

- Provides stable image ✓
- No power or hydraulics ✓
- Easy and flexible ✓
- Free up assistant ✓

